
# Način primene AQUACEL® Ag\*

## Aplikacija\*

- Ranu isprati vodom ili fiziološkim rastvorom
- Postaviti oblogu na površinsku ranu tako da prelazi ivice rane za najmanje 1 – 2 cm (Slika 1)
- Duboku ranu ispuniti sa oblogom ili trakom do 80% dubine zbog bubrenja obloge, a ostatak obloge ili trake ostaviti da viri 2,5 cm iznad rane kako bi se lakše uklonila sa rane (Slika 2)
- Pokriti sa sekundarnim zavojnim materijalom (npr. Granuflex extra thin) ili gazom u zavisnosti od tipa rane

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
### Reference:

<sup>1</sup>Harding KG, Price P, Robinson B, Thomas S, Hofman D. Cost and dressing evaluation of hydrofiber and alginate dressing in the management of community-based patients with chronic leg ulceration. Wounds 2001;13:229-236. <sup>2</sup>Data on File:MA068, ConvaTec. <sup>3</sup>American Diabetes Association. Consensus development conference on diabetic foot wound care. Diabetes Care 1999;22:1354-1360. <sup>4</sup>Data on File:Alt11, ConvaTec. <sup>5</sup>Foster L, Moore P. The application of a cellulose-based fibre dressing in surgical wounds. J Wound Care. 1997;6:469-473. <sup>6</sup>Bowler PG, Jones SA, Davies BJ, Coyle E. Infection control properties of some wound dressings. J Wound Care 1999;8:499-502. <sup>7</sup>Data on File:WA129, ConvaTec. <sup>8</sup>Foster L, Moore P, Clark S. A comparison of hydrofiber and alginate dressings in the management of deep ulcerations of diabetic foot. Diabetes Med. 2001;18:320-324. <sup>10</sup>Chen J, Hutchinson J, Waring M, et al. AQUACEL® Ag Hydrofiber™ wound dressing. ConvaTec/Bristol-Myers Squibb; 1998;7. <sup>11</sup>Robinson BJ. The use of a hydrofiber dressing in wound management. J Wound Care. 2000;9:32-34. <sup>12</sup>Armstrong SH, Ruckley CV. Use of a fibrous dressing in exuding leg ulcers. J Wound Care. 1997;8:322-324. <sup>13</sup>Data on File:MA067, ConvaTec. <sup>14</sup>Demling RH, DeSanti L. Effects of silver on wound management. Wounds 2001;13 (suppl):A5-A15. <sup>15</sup>Data on File:MA072, ConvaTec. <sup>16</sup>Data on File:Alt10, ConvaTec. <sup>17</sup>Sibbald RG, Williamson D, Orsted HL, et al. Preparing the wound bed – debridement, bacterial balance, and moisture balance. Ostomy Wound Manage. 2000;46:14-20,22,24-28, 30-35. <sup>18</sup>Lydson MJ. The development of Aquacel® Hydrofiber™ dressings. In: Krieg T, Harding KG, eds. AQUACEL® Hydrofiber™ dressing: the next step in wound dressing technology. Proceedings of a Satellite Symposium at the 6th Congress of the European Academy of Dermatology & Venerology, Dublin, Ireland, 11-15 September 1997. London, UK: Churchill Communications; 1998:1-3. <sup>19</sup>Data on File:WA 136, ConvaTec.

The Queen's Award  
for Technological Achievement

ConvaTec  
A Bristol-Myers Squibb Company

AQUACEL® Ag  
Hydrofiber™ Dressing with Silver

TT MEDIK

TT Medik, distributer za ConvaTec  
Bulevar Mihajla Pupina 10D/1  
11070 Novi Beograd, Srbija  
Tel: 011 311 51 52, 311 51 53  
office@ttmedik.co.yu  
www.ttmedik.co.yu

Nežan spolja<sup>1</sup>  
Moćan iznutra<sup>2</sup>

AQUACEL® Ag  
Hydrofiber™ Dressing with Silver

Antimikrobna obloga sa srebrom


„Pre upotrebe detaljno proučiti uputstvo! O indikacijama, riziku upotrebe i neželjenim dejstvima na medicinsko sredstvo posavetujte se sa lekarom ili farmaceutom.”

AQUACEL® Ag  
Hydrofiber™ Dressing with Silver

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