

Samo za stručnu javnost

**NOVO!**

# *Foam*Lite™

ConvaTec

ConvaTec program za negu rana  
Adhezivna penasta silikonska obloga



**TT MEDIK**

Nosilac upisa medicinskog sredstva u Registar  
Broj i datum rešenja o upisu medicinskog sredstva u Registar:  
515-02-01468-18-001 od 13.12.2018.

 **ConvaTec**


## FoamLite™ ConvaTec obloga posebno je dizajnirana za tretman akutnih i hroničnih rana bez eksudata ili sa malom količinom eksudata.


FoamLite™ ConvaTec je tanka penasta obloga sa filmom koji omogućava rani da „diše“ i sadrži perforirani nežni silikonski adheziv koji prijta koži. \* Ova obloga je fleksibilna i lako se prilagođava svim delovima tela.


**1 Vodootporna barijera**, mekana na dodir, štiti ranu od bakterija i virusa.†


**2 Tanki i fleksibilni sloj poliuretanske pene** upija malu količinu eksudata i održava vlažnu sredinu pogodnu za zarastanje rane.†1

**3 Perforirani, nežni silikonski adheziv, prijatan za kožu.** Omogućava jednostavnu aplikaciju, repozicioniranje i atraumatsko uklanjanje, bez oštećenja rane i okolne kože.




 Fleksibilan i prilagodljiv†1

 Barijera za bakterije i viruse†1

 Može se seći i oblikovati\*


 Adheziv se ne lepi za rukavice ili za druge delove obloge†

 Maksimalno preporučeno vreme stajanje obloge na rani je 7 dana\*


 Jednostavan za primenu, repozicioniranje i uklanjanje†1,2

 Vodootporan/pogodan za tuširanje†1


## Indikacije za primenu






1. Akutne rane



2. Hronične rane



3. Hirurške rane



† Demonstrirano in vitro

\* Pogledajte uputstvo za upotrebu FoamLite™ ConvaTec obloge za kompletnu informaciju o načinu upotrebe




**FoamLite™**  
ConvaTec



FoamLite™ ConvaTec obloga za tretman rana sa malom količinom eksudata ili bez eksudata

Veličina	Pakovanje	Kat. broj
<b>Adhezivna</b>		
5.5 x 12cm	10	421563
8 x 8cm	10	421557
10 x 10cm	10	421559
15 x 15cm	10	421561



1. *In vitro* Performance of Foam Lite™ ConvaTec WHRI4680 MS132. 25th April 2016. 2. Meuleneire F, Rücknagel H. Soft silicones Made Easy. Wounds International 2013 (May). Available from: [www.woundsinternational.com](http://www.woundsinternational.com)

**TT MEDIK**

 **ConvaTec**

TT Medik, distributer za ConvaTec, Bulevar Mihajla Pupina 10D/1,  
11070 Novi Beograd, Srbija, Tel: 011 311 51 52, 311 51 53  
[office@ttmedik.co.rs](mailto:office@ttmedik.co.rs), [www.ttmedik.co.rs](http://www.ttmedik.co.rs)